

Product Datasheet

Recombinant Mouse Stromal cell-derived factor 2-like protein 1(Sdf2l1) (orb1631310)

Description	Recombinant Mouse Stromal cell-derived factor 2-like protein 1(Sdf2l1)
Species/Host	Yeast
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 6xHis-tagged and C-terminal Myc-tagged
Note	For research use only
Expiration Date	6 months from date of receipt.