


Product Datasheet

ASMTL antibody (orb1634)

Description	ASMTL antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ASMTL (31-130/621 aa)
Target	ASMTL
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	69 kDa
Uniprot ID	O95671
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human colon tissue showing ASMTL expression.