

Product Datasheet

Human PARVA protein (orb163434)

Description	Recombinant human PARVA protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 50mM Tris, 150mM NaCl, 40% Glycerol, pH 7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of 50mM Tris-HCl, 150mM NaCl, 40% Glycerol, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Met1-Glu372
Source	E.coli
Uniprot ID	Q9NVD7
Expiration Date	6 months from date of receipt.