

Product Datasheet

Human DFFA protein (orb163435)

Description	Recombinant human DFFA protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Protein Sequence	Met1-Thr331
Source	E.coli
Uniprot ID	AAH07721.1
Expiration Date	6 months from date of receipt.