

Product Datasheet

Human METAP1D protein (orb163445)

Description	Recombinant human METAP1D protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 50mM Tris, 100mM NaCl, pH 8.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mMTris-HCl,pH8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His&C-His
Note	For research use only
Application notes	6His tag at the N-terminus, 6His tag at the C-terminus.
Protein Sequence	Arg44-Ala335
Source	E.coli
Uniprot ID	Q6UB28
Expiration Date	6 months from date of receipt.