

Product Datasheet

Human DCK protein (orb163448)

Description	Recombinant human DCK protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl pH7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His&T7
Note	For research use only
Application notes	6His, T7 tag at the N-terminus.
Protein Sequence	Met1-Leu260
Source	E.coli
Uniprot ID	P27707
Expiration Date	6 months from date of receipt.