

Product Datasheet

Human RPE protein (orb163474)

Description	Recombinant human RPE protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 6.2.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 6.2.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Met1-Arg228
Source	E.coli
Uniprot ID	Q96AT9
Expiration Date	6 months from date of receipt.