

Product Datasheet

Human BPGM protein (orb163476)

Description	Recombinant human BPGM protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris, 1mM DTT, pH 8.0 .
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 1mM DTT, pH 8.0 .
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Ser2-Lys259
Source	E.coli
Uniprot ID	P07738
Expiration Date	6 months from date of receipt.