

Product Datasheet

Human PHLDA2 protein (orb163482)

Description	Recombinant human PHLDA2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 0.1M NaCl, 1mM DTT, pH 8.0.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mMTris-HCl,350mMNaCl,1mM DTT,pH8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Met1-Pro152
Source	E.coli
Uniprot ID	Q53GA4
Expiration Date	6 months from date of receipt.