

Product Datasheet

Human ACADM protein (orb163483)

Description	Recombinant human ACADM protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris, 0.1M NaCl, 20% Glycerol, pH 8.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 100mM NaCl, 20% Glycerol, pH 8.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Protein Sequence	Lys26-Asn421
Source	E.coli
Uniprot ID	P11310
Expiration Date	6 months from date of receipt.