

Product Datasheet

Human Blk protein (orb163491)

Description	Recombinant human Blk protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris, 500mM NaCl, 1mM DTT, pH 7.4.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 500mM NaCl, 1mM DTT, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Gly2-Pro505
Source	E.coli
Uniprot ID	P51451
Expiration Date	6 months from date of receipt.