

Product Datasheet

Human CAMK2B protein (orb163492)

Description	Recombinant human CAMK2B protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 10mM HEPES, 150mM NaCl, pH 7.5.
Form/Appearance	Supplied as a 0.2 µm filtered solution of 10mM HEPES, 150mM NaCl, pH 7.5.
Storage	Store at ≤-70°C, stable for 6 months after receipt.Store at ≤-70°C, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Met1-Gln503
Source	E. coli
Uniprot ID	Q13554
Expiration Date	6 months from date of receipt.