

Product Datasheet

Human Glutathione S-Transferase protein (orb163494)

Description	Recombinant human Glutathione S-Transferase protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 20mM Tris, 200mM NaCl, pH 7.0.
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 200mM NaCl, pH 7.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Protein Sequence	Met1-Leu199
Source	E.coli
Uniprot ID	O60760
Expiration Date	6 months from date of receipt.