

Product Datasheet

Human Caspase-10 protein (orb163519)

Description	Recombinant human Caspase-10 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Supplied as a 0.2 µm filtered solution of 25mM HEPES, 10mM DTT, pH 7.5.
Form/Appearance	Supplied as a 0.2 um filtered solution of 25mM HEPES, 10mM DTT, pH 7.5.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Val220-Ile480
Source	E.coli
Uniprot ID	Q92851
Expiration Date	6 months from date of receipt.