

## Product Datasheet

### Human SULT1A3 protein (orb163526)

|                          |                                                                                        |
|--------------------------|----------------------------------------------------------------------------------------|
| <b>Description</b>       | Recombinant human SULT1A3 protein                                                      |
| <b>Endotoxins</b>        | Less than 0.1 ng/ug (1 EU/ug)                                                          |
| <b>Preservatives</b>     | Supplied as a 0.2 µm filtered solution of 20mM Tris,100mM NaCl, pH 8.0.                |
| <b>Form/Appearance</b>   | Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 100mM NaCl, pH 8.0.           |
| <b>Storage</b>           | Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles. |
| <b>Tag</b>               | N-His                                                                                  |
| <b>Note</b>              | For research use only                                                                  |
| <b>Application notes</b> | 6His tag at the N-terminus.                                                            |
| <b>Protein Sequence</b>  | Met1-Leu295                                                                            |
| <b>Source</b>            | E.coli                                                                                 |
| <b>Uniprot ID</b>        | <b>P50224</b>                                                                          |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                         |