

Product Datasheet

SARS-CoV-2 (2019-nCoV) Spike Protein (S1 Subunit, mFc Tag) (orb1635308)

Description	SARS-CoV-2 (2019-nCoV) Spike Protein (S1 Subunit, mFc Tag)
Species/Host	Human
Storage	Samples are stable for up to twelve months from date of receipt at -20°C to -80°C Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.