

Product Datasheet

Human Myozenin-2 protein (orb163531)

Description	Recombinant human Myozenin-2 protein
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Preservatives	Lyophilized from a 0.2 µm filtered solution of 10mM Tris, pH 8.0.
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 10mM Tris-HCl, pH 8.0.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Protein Sequence	Met1-Leu264
Source	E.coli
Uniprot ID	Q9NPC6
Expiration Date	6 months from date of receipt.