

Product Datasheet

SOX8 antibody (orb163568)

Description	Rabbit polyclonal antibody to SOX8
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SOX8
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47314
Uniprot ID	P57073
NCBI	NP_055402.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.