

Product Datasheet

TBX15 antibody (orb163578)

Description	Rabbit polyclonal antibody to TBX15 (N-term).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	TBX15
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	65757
Uniprot ID	Q96SF7
NCBI	NP_689593.2
Dilution Range	IF/ICC: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.