

Product Datasheet

HER2 (phospho-Y1127) antibody (orb163601)

Description	Rabbit polyclonal antibody to HER2 (phospho-Y1127)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT, IF, WB
Immunogen	Synthetic Peptide
Target	ERBB2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	137910
Uniprot ID	P04626
NCBI	NP_004439.2 , NP_001005862.1
Dilution Range	WB: 1:5000, IF/ICC: 1:100, DOT: 1:500
Expiration Date	12 months from date of receipt.