

Product Datasheet

AMPD3 antibody (orb163607)

Description	Rabbit polyclonal antibody to AMPD3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	AMPD3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	88812
Uniprot ID	Q01432
NCBI	NP_000471.1 , NP_001165901.1 , NP_001020560.1 , NP_001165902.1 , NP_001020561.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.