

Product Datasheet

ASAH2 antibody (orb163661)

Description	Rabbit polyclonal antibody to ASAH2 (C-term).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ASAH2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	85516
Uniprot ID	Q9NR71
NCBI	NP_001137446.1 , NP_063946.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.