

Product Datasheet

SNAP25 (phospho-T138) antibody (orb163730)

Description	Rabbit polyclonal antibody to SNAP254 which is also known as SUP whose expression can be detected in...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	DOT
Immunogen	Synthetic Peptide
Target	SNAP25
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23315
Uniprot ID	P60880
NCBI	NP_570824.1, NP_003072.2
Dilution Range	DOT: 1:500
Expiration Date	12 months from date of receipt.