

Product Datasheet

S. Aureus Delta Hemolysin antibody (orb163731)

Description	Rabbit polyclonal antibody to S. Aureus Delta Hemolysin.
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	hld
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	2979
Uniprot ID	P0C1V1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

