




Product Datasheet

COX7A1 antibody (orb163736)

Description	Rabbit polyclonal antibody to COX7A1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	COX7A1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	9118
Uniprot ID	P24310
NCBI	NP_001855.1
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of human skeletal ...



Western blot analysis of human skeletal ...