

Product Datasheet

KANK2 antibody (orb163752)

Description	Rabbit polyclonal antibody to KANK2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KANK2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	91174
Uniprot ID	Q63ZY3
NCBI	NP_001129663.1 , NP_056308.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.