

Product Datasheet

NEIL2 antibody (orb163769)

Description	Rabbit polyclonal antibody to NEIL2 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	NEIL2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36826
Uniprot ID	Q969S2
NCBI	NP_001129220.1 , NP_001129219.1 , NP_001129218.1 , NP_659480.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.