

Product Datasheet

RPS14 antibody (orb163771)

Description	Rabbit polyclonal antibody to RPS14 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	RPS14
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	16273
Uniprot ID	P62263
NCBI	NP_005608.1, NP_001020241.1, NP_001020242.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.