

Product Datasheet

P2RY6 antibody (orb163778)

Description	Rabbit polyclonal antibody to P2RY6 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	P2RY6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	36429
Uniprot ID	Q15077
NCBI	NP_001264137.1 , NP_001264136.1 , NP_001264134.1 , NP_001264133.1 , NP_789767.1 , NP_789768.1 , NP_789766.1 , NP_001264135.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.