

Product Datasheet

SNX17 antibody (orb163788)

Description	Rabbit polyclonal antibody to SNX17 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SNX17
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52901
Uniprot ID	Q15036
NCBI	NP_001253989.1 , NP_055563.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.