

Product Datasheet

Brsk2 antibody (orb163790)

Description	Rabbit polyclonal antibody to Brsk2
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	Brsk2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	81733
Uniprot ID	Q69Z98
NCBI	NP_001009929.1 , NP_001009930.1 , NP_083702.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.