

Product Datasheet

PMP22 antibody (orb163801)

Description	Rabbit polyclonal antibody to PMP22.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PMP22
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	17891
Uniprot ID	Q01453
NCBI	NP_000295.1 , NP_696996.1 , NP_001268384.1 , NP_696997.1 , NP_001268385.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

