

Product Datasheet

GABPB1 antibody (orb163835)

Description	Rabbit polyclonal antibody to GABPB1 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GABPB1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42483
Uniprot ID	Q06547
NCBI	NP_057738.1 , NP_002032.2 , NP_005245.2 , NP_057739.1 , NP_852092.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.