

Product Datasheet

NR2C1 antibody (orb163850)

Description	Rabbit polyclonal antibody to NR2C1 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	NR2C1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	67315
Uniprot ID	P13056
NCBI	NP_001027458.1 , NP_001120834.1 , NP_003288.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.