

Product Datasheet

TTC23 antibody (orb163868)

Description	Rabbit polyclonal antibody to TTC23
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TTC23
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	50009
Uniprot ID	Q5W5X9
NCBI	NP_001035750.1 , NP_001035748.1 , NP_001035749.1 , NP_001035746.1 , NP_001035745.1 , NP_001035747.1 , NP_075056.3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.