

Product Datasheet

NCOA4 antibody (orb163902)

Description	Rabbit polyclonal antibody to NCOA4 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	NCOA4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	69726
Uniprot ID	Q13772
NCBI	NP_001138734.1 , NP_001138732.1 , NP_001138735.1 , NP_005428.1 , NP_001138733.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.