

Product Datasheet

L1CAM antibody (orb163906)

Description	Rabbit polyclonal antibody to L1CAM (C-term).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	L1CAM
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	140003
Uniprot ID	P32004
NCBI	NP_076493.1, NP_001137435.1, NP_000416.1, NP_001265045.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.