

Product Datasheet

ESR1 antibody (orb163907)

Description	Rabbit polyclonal antibody to ESR1.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	ESR1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	66216
Uniprot ID	P03372
NCBI	NP_001116212.1 , NP_001116213.1 , NP_000116.2 , NP_001116214.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

