

Product Datasheet

CBFA2T2 antibody (orb163927)

Description	Rabbit polyclonal antibody to CBFA2T2 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	CBFA2T2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	67133
Uniprot ID	O43439
NCBI	NP_001028171.1 , NP_001034798.1 , NP_005084.1
Dilution Range	IF/ICC: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.