

Product Datasheet

ZNF512 antibody (orb164031)

Description	Rabbit polyclonal antibody to ZNF512 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ZNF512
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	64682
Uniprot ID	Q96ME7
NCBI	NP_001258216.1 , NP_001258218.1 , NP_001258215.1 , NP_001258247.1 , NP_115810.2 , NP_001258217.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.