

Product Datasheet

HTR7 antibody (orb164037)

Description	Rabbit polyclonal antibody to HTR7 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	HTR7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53555
Uniprot ID	P34969
NCBI	NP_062874.1 , NP_062873.1 , NP_000863.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.