

Product Datasheet

CRIP2 antibody (orb164114)

Description	Rabbit polyclonal antibody to CRIP2 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CRIP2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22493
Uniprot ID	P52943
NCBI	NP_001303.1, NP_001257766.1, NP_001257770.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.