

Product Datasheet

ZC3H14 antibody (orb164128)

Description	Rabbit polyclonal antibody to ZC3H14
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ZC3H14
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	82876
Uniprot ID	Q6PJT7
NCBI	NP_997545.2 , NP_997543.1 , NP_001153576.1 , NP_997544.1 , NP_001153575.1 , NP_079100.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.