

Product Datasheet

NKIRAS2 antibody (orb164136)

Description	Rabbit polyclonal antibody to NKIRAS2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	NKIRAS2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	21508
Uniprot ID	Q9NYR9
NCBI	NP_001138399.1 , NP_060065.2 , NP_001138400.1 , NP_001001349.1 , NP_001138401.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.