

Product Datasheet

SH3PXD2A antibody (orb164168)

Description	Rabbit polyclonal antibody to SH3PXD2A (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SH3PXD2A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	125289
Uniprot ID	Q5TCZ1
NCBI	NP_055446.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.