

Product Datasheet

FAIM antibody (orb164212)

Description	Rabbit polyclonal antibody to FAIM
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	FAIM
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	20215
Uniprot ID	Q9NVQ4
NCBI	NP_001028204.1 , NP_001028203.1 , NP_001028202.1 , NP_060617.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.