

Product Datasheet

CPEB1 antibody (orb164216)

Description	Rabbit polyclonal antibody to CPEB1 (C-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CPEB1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	62595
Uniprot ID	Q9BZB8
NCBI	NP_001073001.1 , NP_001073002.1 , NP_085097.3 , NP_001073003.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.