

Product Datasheet

DTNBP1 antibody (orb164252)

Description	Rabbit polyclonal antibody to DTNBP1 (N-term).
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	DTNBP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	39493
Uniprot ID	Q96EV8
NCBI	NP_115498.2 , NP_001258598.1 , NP_001258597.1 , NP_898861.1 , NP_001258596.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.