

Product Datasheet

SH3BP2 antibody (orb164254)

Description	Rabbit polyclonal antibody to SH3BP2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SH3BP2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	62244
Uniprot ID	P78314
NCBI	NP_003014.3, NP_001139327.1, NP_001116153.1, NP_001139328.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.