

Product Datasheet

CD79B antibody (orb164260)

Description	Rabbit polyclonal antibody to CD79B.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CD79B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	26048
Uniprot ID	P40259
NCBI	NP_001035022.1 , NP_067613.1 , NP_000617.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.